

**ORDER YOURS
TODAY AT
SheRexRecovery.com**

Questions?

Email: info@SheRexRecovery.com

Call: 303-356-7249





HEAL
TODAY
IMAGINE
TOMORROW



SHE REX[®]
recover strong.


WELCOME!! YOU'VE JUST HAD BREAST SURGERY...

...and you're looking forward to getting back to all of life's activities. No matter why you decided to go on this journey, the next few weeks will include some emotional and physical changes. Although people arrive at the decision to undergo breast surgery from many different paths, something everyone who completes breast augmentation or reduction can agree upon is a desire for a quicker recovery with less setbacks to healing.

It's time to commit to a smooth recovery. Give yourself permission to more easily take care of yourself while going about your daily activities and to focus on the things most important to you, without worrying about accidentally overextending a motion. That's what **SHE REX®** is all about, seamlessly and comfortably integrating into your routine – gently reminding you to keep movements within the best range for optimal recovery.

Created by Kelli Beckish, a successful entrepreneur and textile manufacturer, after her own breast surgery, SHE REX® is designed to make the healing process more efficient. The product is an extraordinarily simple, but truly ingenious invention to help those recovering keep their movements limited to the safe range. The name **SHE REX®** is an acknowledgment of how the limited movements make the user look a bit like the mighty dinosaur with short arms. Today, **SHE REX®** is available to help anyone recovering from augmentation, reduction, or other breast surgery procedures.





GENTLE RECOVERY REMINDER

SHE REX® patent-pending post-op wearable bands comfortably remind you to keep your arms from overreaching – avoiding the painful movements that hinder recovery. Ideal for use following mastectomies, augmentations and other surgical breast procedures.

HOW TO ADJUST SHE REX®



A

B*

C

Adjust buckle **A** to fit right or left arm.

Adjust buckle **B*** to size around your back.

Adjust buckle **C** to fit opposite arm.

Band sizes are comparable to T-shirt sizes. There's 6 extra inches of strap for a proper adjustable fit. * buckle B is available in sizes L-2XL

